# Evaluation of the potential to cause dermal irritation of the Samples provided by the Customer

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## Utilizing the *Utilizzando il* Dermal Irritection® test method

### **Enclosed (*Allegati*):**

- Introduction
- Scientific Background
- Bibliography
- Information for the customized service

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### **the reasons why using the IRRITECTION™ Assay System.**

The Irritection™ Assay System (IAS) consists of newly-developed test kits, instrumentation and computer software that have been integrated to provide an automated in vitro testing capability that detects, predicts and ranks the ocular and Dermal irritation potential of several types of materials. Because the system is reliable, rapid and cost-effective, it can be readily employed as a broad-screen complement or as a complete replacement test for in vivo test methods. Consequently, the Irritection system can be used as an efficient tool that substantially reduces the cost of developing products such as cosmetics, pharmaceuticals, petrochemical formulations, surfactant-containing compounds and a wide variety of consumer products. In the proposed OECD testing and evaluating strategy for eye irritation/corrosion (ECVAM Workshop Report 34, page 64) Ocular Irritection™ could be the single valid test which encloses the assessments of point 5 (severe eye irritant) and point 6 (eye irritation). In case of validation, only assessment of skin corrosion (point 7) would be necessary, with no further tests on rabbit eye.

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INT.E.G.R.A. S.r.l. Introduction

## SCIENTIFIC BACKGROUND

The Irritection® test system is an innovative test that includes the use of purposely designed matrix and of the equipment necessary to read and elaborate the resulting data with a dedicated software that permits the identification and prediction of the dermal irritation potential of the products in the different forms available.


This system is extremely fast and reliable and therefore can be used to add the irritant factor to previously completed product experiments or to extend the tests made on larger ranges of products.

The initial reason for developing the Irritection® test was the idea that an *in vitro* test could reliably substitute the Patch test. Consequently, the research focused on the understanding of the mechanisms that are responsible for the corneous opacification and the skin irritation observed in *in vivo* tests with the intention of identifying the matrix that could substitute the dermal and corneous protein damaged through contact with the substances under examination. This research brought about the development of the Eyetex™ and Skintex™ protocols that in turn led to the development of the current Irritection® test.

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**Figure 1.** The Dermal Irritection Model (*Modello operativo del protocollo*)

It is a known fact that eye and skin irritation tests are based on well known biochemical mechanisms. In particular, it has been demonstrated that the corneous, in normal

DI001

INT.E.G.RA. S.r.l. Introduction

conditions, is transparent due to structural formation of the proteins and the carbohydrates form which it is formed and which show a marked state of internal organisation. When the irritating agents are placed on the eye, they modify the quaternary structural formation of the proteins, causing the corneous opacification that constitutes a parameter of observation of the Draize test. In the development of the Irritection® test, it was discovered that this process can be emulated *in vitro* when the irritant agents are applied to a mixture of proteins and glycoproteins.

The main component of the reagent matrix used in the Irritection® test is an oligomeric protein formed of 12 sub-units, with an apparent molecular weight of around 320.000 D. The proteic oligomers tend to form an auto-associative structure, forming associations bonded by non-covalent interactions. Such a matrix also contains carbohydrates, lipids and components with a low molecular weight, which, in mixtures with other components in a physiological solution, reproduce a similar structure to that of the skin.

It has been demonstrated that when irritant agents are applied to the matrix of the Irritection® test, they produce very similar effects to those occurring *in vivo*. Furthermore, the effects of interference upon the proteic structure that interfere both with the corneous *in vivo* and with the artificial matrix of the Irritection® test can be correlated to the aggressiveness of the irritant agent, in that the nucleation of particles that make the matrix cloudy result to be proportionate to the irritant factor of the tested substance. This phenomenon can be quantified with a optical density spectrophotometer measurement performed with a custom prepared spectrophotometer. According to the Dermal Irritection® protocol, the optical density measurements are performed at 450nm (OD<sub>450</sub>), while according to the Ocular Irritection® protocol they are performed at 400nm (OD<sub>400</sub>).

Due to the fact that the length of time that the substances under examination are in contact with the matrix is a critical parameter, a semi-permeable membrane with controlled porosity is used which makes the identification of the maximum opacification point much easier. The point of maximum opacification is then determined by spectrophotometer.

Based on the tests performed in parallel to determine the calibration curve, it was possible to identify a direct correlation between the Patch test and the scale determined with the Irritection® protocol (Irritection Score), according to the following summary:

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Tab.1 Relationship of Irritection Equivalent Score to Irritancy Classification for the Dermal Irritection® Test Method.  
Corrispondenza fra scala Dermal Irritection® e Irritazione cutanea prevista

DI001

INT.E.G.R.A. S.r.l. Introduction

### **Interpretation of results**

The optical density values elaborated by the software at the interface with a UV spectrophotometer are divided into 3 groups: **Sample values, Calibrator Values e Quality Control Values**. The Quality Control Values are formed of two samples contained within each Irritection® test kit, and which supply an irritant score in the two main windows (irritant and non irritant), with a area of tolerance. The substances have been chosen from amongst those that cause the greatest difficulty for the test so that a self-test is created. If any of the parameters results to be outside the foreseen tolerances, the test results as non-qualified, and the software does not give the elaborated results.

The Calibrators are 4 samples that are supplied with each kit and have the double function of validating the batch of the matrix/porous septum used, and of supplying the best correlation line for converting the Optical Density (OD) into an Irritant Score (IDE) for the batch of the kit used. The 4 Calibrator samples are naturally in the 4 windows, and if one of the samples does not give a result within the programmed tolerance, the test results as not-qualified and the software does not give the elaborated results.

According to the OD obtained with the Calibrators, if the Quality Controls are within tolerance the software correlates the OD of the samples with the equivalent irritant score of the Patch test. The most significant value, based on the correlated tests performed in parallel with the Patch test, proves to be the maximum value amongst those evaluated (Maximum Qualified Score).

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INT.E.G.R.A. S.r.l. Introduction

## Bibliography

1. Draize, J.H., Woodard, G., Calvery, H.O. Methods for the study of irritation and toxicity of substances applied topically to the skin and mucous membranes. *J. Pharmacol. Exp. Ther.* 82, 377-390 (1944).
2. Kelly, C. EYETEX : An in-vitro method of predicting ocular safety. *Pharmacopeial Forum*, 4815-4824 (1989).
3. Gordon, V.C. The scientific basis of the EYETEX<sup>TM</sup> system. *ATLA* 20, 537-548 (1992).
4. Cooper, J.A. Saracci, R. and Cole, P. Describing the validity of carcinogen screening tests. *Br.J. Cancer*, 39, 87-89 (1979).
5. Scala, R.A. Theoretical approaches to validation. In : Alternative Methods in toxicology, vol.5, *In Vitro toxicology: Mechanisms and New Technology* (Ed. A.M. Goldberg) pp. 1-9, Mary Ann Liebert Publishing, New York, New York (1987).
6. Roy, T.A., Saladdin, K.A., Yang, J.J. and Gordon, V.C. Evaluation of the EYETEX<sup>TM</sup> system as screen for the eye irritancy of petroleum products. In: Alternative Methods in Toxicology, vol.8, *In Vitro Toxicology: Mechanisms and New Technology* (Ed. A.M. Goldberg) pp. 223-231, Mary Ann Liebert Publishing, New York, New York (1991).
7. Kruszewski, F.H., Hearn, L.H., Smith, K.T., Teal, J.J., Gordon, V.C. and Dickens, M.S. Application of the EYETEX<sup>TM</sup> System to the evaluation of cosmetic products and their ingredients. *ATLA*, 20, 146-163 (1992).
8. Sina, J.F., Galer, D.M., Sussman, R.G., Gautheron, P.D., Sargent, E.V., Leong, B., Shah, P.V., Curren, R.D. and Miller, K. A collaborative evaluation of seven alternatives to the Draize eye irritation test using pharmaceutical intermediates. *Fundamental and Appl. Toxicology*, 26, 20-31 (1995).
9. Règnier, J-F., Imbert, C. and Boutonnet, J-C. Evaluation of the EYETEX<sup>TM</sup> System as a screening for the ocular irritancy of cosmetic product. *ATLA*, 22, 32-50 (1994).

DI001

INT.E.G.R.A. S.r.l. Introduction

To INT.E.G.R.A. CUSTOMIZED SERVICE

**Required Information concerning the sample:**

	Sample N.:	Sample N.:	Sample N.:	Sample N.:
Sample/s name (code, batch):				
Product/s Classification (Surfactant?):				
Solubility:				
Waste Class (when applicable):				
Storage information:				
Please mark the following information for each product, when it is "true", with X				
<i>Acidic Formulations pH &lt; 2.0</i>				
<i>Alkaline Formulations pH &gt; 9.0</i>				
<i>Intensely colored materials</i>				
<i>Oily and water-insoluble organic chemicals</i>				
<i>Volatile Ketones &gt; 10%</i>				
<i>Free Nonionic Surfactants &gt; 10%</i>				
<i>Sorbitol at concentrations &gt; 5%</i>				
<i>Urea at concentrations &gt; 5%</i>				
<i>Manganese Violet</i>				
<i>Aluminium Chlorohydrate</i>				
<i>Aluminium Zirconium Chlorohydrate</i>				
<i>Aluminium Chloride</i>				
<i>Titanium oxide &gt; 10%</i>				
<i>Zinc Oxide &gt; 10%</i>				
<i>Silver Salts contained</i>				
<i>Ferrous sulfate contained</i>				
<i>Zinc Sulfate contained</i>				

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[WWW.cosmeting.it](http://WWW.cosmeting.it)

(please send this information by fax or e-mail)

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INT.E.G.R.A. S.r.l.: 20056 Via G. Pastore, 3 Trezzo Sull'Adda (MI) Italy